

## User's Manual

Description for VSZT-1 bioetheric field harmonizer (Orgon-radiator) An efficient tool for natural treatments at home.

VSZT-1 device is a special energetic and bioinformatics system. The device is an environment-powered unit that extracts energy for its operation from the surrounding space. Therefore it is very cost efficient. Its active parts are charged crystals, and the chassis is made of a special aluminum alloy. The device does not contain any moving parts of electronics; therefore the chance of a malfunction is very minimal. The operating life of a VSZT-1 can reach or exceed 100 years unless the tool is damaged or gets dismantled. There is an opportunity to encode the crystals with special programs (e.g. for veterinary use).

**Usage of VSZT-1 does not replace medical treatments and does not entitle the user to stop taking prescribed medications or to alter the medical therapy. Please consider carefully! Continue medical visits as told!**

Most important areas of application of VSZT-1:

- 1) Supplementary therapy for tumor treatments.
- 2) Rehabilitation after cardiac arrests and other cardio-vascular problems.
- 3) Rehabilitation after brain damage. Children born with low oxygen shock.
- 4) Improving general conditions of patients after surgeries. Enhancing the healing processes of wounds and preventing side-effects.
- 5) Treating inhaling problems, allergies.
- 6) Treating infertility, endangered childbirth, gynecological problems.
- 7) Athlete's foot, vein related problems.
- 8) Refreshing the body in exhausted or immuno-deficiency states.
- 9) Enhancing the digestive systems.
- 10) Decreasing the effects of climax.
- 11) Joint related complaints.
- 12) Epilepsy, phobias, headache.
- 13) SM, Parkinson's disease.
- 14) Prostate related illnesses, urination problems.
- 15) Hair related problems.
- 16) Against osteoporosis.
- 17) Supplementary treatment for diabetes.
- 18) Supplementary treatment for hypotension.

**In its standard setup the device is not suitable for treating mental illnesses!**

For use in energy deficiency state caused by earth radiation, for the replacement of energy, and for the provision of the defense against harmful earth radiation.

## Installation instructions

### Parts:

- Radiator head
- Information jar
- Head-locking screw
- Cable between head and tube
- Stander screw
- Stand
- Base

### Contents:

- 1 carrying case
- 1 VSZT device in parts
- 1 User's guide

**In its default setting the lock is blocked, which can be activated after the removal of the small plastic piece from the case's main lock.**

### ***Assembling instructions:***

We start with the stander. Take the base out from the carrying case. Put the screw into the hole and insert it into the stand's borehole. Drive it in and fasten the stander as is. Unscrew the jagged head on top of the stander and insert the rocket shaped main radiator's side pin into the hole on the side of the column. Fasten the screw and adjust angle of head. Insert one end of the cable into the bottom side of the information jar and the other end into the back of the radiator head.

If not in use, VSZT-1 should be kept in the carrying case or in a charging position. **Charging is done in a way that the information jar is placed in front of the radiator head so that the tip of the head touches the brass colored top of the information jar, thus it creates a short circuit with the cable.** It is advisable to keep the device in such a position for at least 20 minutes a day. Charging time may exceed even days or months without any damage done to the device.

**Beware of the sharp tip of the radiator head during treatments so that it does not render any damage to living tissues. Keep it locked away from kids. Only apply to kids with parental attention.**

**Clothing does not deter the radiation of the device, can be applied without the need to remove clothes.**

### Operation:

- For general treatment, use the device for 30-60 minutes.
- For traumatic treatment or severe pains, apply it for 15-30 minutes directly to the target area.

### Suggested usage:

- Prevention of diseases.
- Treating diseases as a supplementary treatment or main one.
- Post-traumatic rehabilitation programs.

### Prevention of diseases:

It can provide extremely efficient prevention especially when there are multiple occurrences of serious hereditary diseases, or a predisposition to these might pose an increased threat to health. Constant energy supplementation ensures good general condition, life processes are able to function harmoniously.

The treatment strengthens the immune system and improves the resistance ability of the body against diseases.

### Use of the appliance in case of different health problems:

In case of diseases related to the head: e.g. / respiratory allergy, headache, cavity inflammations, ear diseases, tingle in the ears, ear inflammations, brain catastrophes, learning difficulties, side effects of oxygen-deficiencies in infants, tumorous diseases, after brain operations, etc.

In case of neck and chest diseases: tonsillitis, swelling of the neck lymph glands, asthma, allergy, breathing difficulties, bronchitis, arhythmia regeneration after stroke, tumours.

In case of abdominal diseases: digestion problems, stomach, gall bladder, intestinal problems, urinary problems, gynaecological, prosthate, abdominal complaints, tumours.

In case of locomotion disorders: strains, fractures, abrasions, inflammations, osteoporosis, rehabilitation.

In case of vascular diseases: narrow blood vessels, vasodilation, discolouration of the capillaries, ulcers, high or low blood pressure, stroke, arhythmia.

### Treatment positions:

In case of the head zone: In a distance of 30 cms from the head, directed at the problematic area Duration of treatment: 15-30 minutes.

In case of the neck and chest: from the adequate direction, from a distance of 20-30 cms, in sitting or lying position. Duration of treatment: 10-60 minutes.

In case of abdominal organs: from the adequate direction, from a distance of 20-30 cms, in sitting or lying position. In case of abdominal complaints e.g. / gynaecological problems, prostate, haemorrhoid problems, menstruation complaints, infertility, etc. / The patient shall lie on her back with her legs slightly apart and the appliance shall be placed between her legs at a distance of 10-15 cms from the body; duration of treatment: 25-30 minutes.

In case of the locomotive organs: the appliance shall be directed at the area to be treated, from a distance of 5-30 cms, for 10-30 minutes.

In case of the blood vessels: the appliance shall be directed at the area to be treated, from a distance of 5-30 cms, for 10-30 minutes per area.

### **Ideas for treating miscellaneous problems**

For prevention: the general treatment target area in this case is the pit of the stomach from a distance of 30-40 cms for 15-30 minutes. Can be repeated daily.

*The treatment of the vegetative nervous system shall be executed here as well.*

*In the case of spinal pain: the painful area is to be exposed directly to the energy radiation from a distance of 8 to 12 cms; duration of treatment: 20-25 minutes. The patient should lie on their side during treatment. The same position is required for treating osteoporosis as well. Treatment of other joints and osteoporosis should observe the same parameters (distance, duration) described above.*

In the case of headaches: Treatment should be directed at the head with a distance of 20-25 cms between the head and the energy transmitter; duration of treatment: 15-30 minutes. Treatment of nervous weakness and learning difficulties should observe the same parameters (distance, duration) described above.

In the case of gynaecological complaints, disorders of menstruation, inflammations, myomas, infertility etc.: The patient must lie on her back with her legs slightly apart and the appliance must be placed between her legs at a distance of 10-15 cms from the body; duration of treatment: 25-30 minutes.

Urination difficulties: frequent urination at night, prostate inflammation, low sperm count, spasmodic states. Haemorrhoids should observe the same parameters (distance, duration) described above.

In the case of toothache: treatment from a distance of 5-10 cms is continued until the pain ceases. Treatment of pain in other parts of the body should observe the same parameters (distance, duration) described above.

Treatment of crural ulcer: The head of the appliance should be pointed at the area to be treated; from a distance of 25-30 cms. The duration of the treatment is at least 30 minutes. The treatment can be repeated several times a day, mainly in case of pains. According to experiences the appliance eases pain effectively. The sensations during treatment, like pins and needles, heat, prickling, contractions are natural phenomena. These sensations and their intensity are not the same for everyone. The ulcer can be treated through the bandage as well, it shall not be removed during treatment.

It can be used in every case when the process of regeneration shall be started. In case of the treatments the treatment plan given by the professional shall be observed, if you do not have such a plan, please choose the above mentioned treatment options.

The distance of the treatment depends on the size of the organ; the bigger the organ, the more the treatment distance shall be. This distance might range between 1 and 30 cms. The duration of treatment is max. 60 minutes. You will not be over charged even if you exceed the treatment period.

For example you fall asleep during treatment, which often happens.

If we place the organ appliance out of this distance, it executes spatial charge and earth radiation protection.

***Please introduce the application of the device gradually. Start with 15 minutes, and increase the treatment time by 5 minutes each week until the maximum duration of treatment as described above. Patients with high blood pressure shall carefully observe the instructions!***

**Warning: if the pain is of an uncertain origin, a doctor must be consulted as the energy treatment suppresses pain which may signal some very serious problem.**

Please use the device creatively. In case of serious diseases it is important to observe the treatment plan developed by a professional. Should you have any questions, please contact the representative of the company or directly to the producer.

**In case of need it is possible to order a jar suitable for the treatment of pets.**

Warranty and after-warranty repairs are done by the service at the following address:  
7025 Bölcské. Gátórház. u. 6. Phone: (+36) 75-535-126 Fax: (+36) 75-535-125  
mobile: (+36) 20 9-391-273

Detailed material regarding the appliance and other useful information on the internet:

Web: [www.kavalko.hu](http://www.kavalko.hu) Web: [www.orgon.hu](http://www.orgon.hu) Mail: [kavalko@t-online.hu](mailto:kavalko@t-online.hu)

#### Warranty

Start of the warranty: 2007 year .....month .....day

Production number: .....

Licence number of the Medical and Hospital Technical Institute: ( ORKI licence number) 541/2000

EU Registration number: HU/CA01/0121/05

Stamp

As a consequence of the special tuning procedure and hygiene reasons we do not redeem any device! The costs of repairs arising from the damage of the device shall be paid by the buyer. Any other faults have a warranty of 1 year. Please notice the above.